

EZ-Retriever[®] System

Dewax, rehydration and antigen retrieval of up to **96** slides in **35** minutes

From the inventors of Antigen Retrieval*



Real Fast
Real Quality
Real EZ



BioGenex

Time Comparison for Deparaffinization + Rehydration + Antigen Retrieval

Conventional System

75 Mins.

Deparaffinization & Rehydration
30 Mins.

Heat Antigen Retrieval/Enzyme
45 Mins.

EZ-Retriever® System

35 Mins.

Deparaffinization, Rehydration &
Antigen Retrieval

Features & Benefits

- Pre-optimized factory protocols
- User-defined protocols
- Precise time and temperature control
- Customized eco-friendly solutions
- Consistent reliable results


Ordering information

EZ-Retriever®	Catalog No.
EZ-Retriever® (110V)	MW014-MO
EZ-Retriever® (230V)	MW013-MO


EZ-AR Elegance Solutions	Catalog No.
EZ-AR™ 1 Elegance	HK546-XAK
EZ-AR™ 2 Elegance	HK547-XAK

EZ-AR™ Solutions	Catalog No.
EZ-AR™ 1	HK521-XAK
EZ-AR™ 2	HK522-XAK
EZ-AR™ 3 (10X)	HK543-YOK
EZ-AR™ 4 (10X)	HK544-YOK
EZ-AR™ Common (5X)	HK545-XOK


Superior Staining




Desmin (AM072) with EZ-AR1 on
Leiomyoma Tissue




Ki-67 (AN527) with EZ-AR2 on
Breast Carcinoma Tissue




Cyclin D1 (AR474) with EZ-AR1 on
Breast Carcinoma Tissue



CDX-2 (AM392) with EZ-AR3 on
Colon Tissue



CD34 (AM236) with EZ-AR1 on
Adeno Carcinoma Tissue



CD30 (AM327) with EZ-AR4 on
Hodgkin Tissue

*Covered by U.S. Patent No. 5,244,787, **Reference:** Shi S-R, Key ME, Kalra KL. Antigen retrieval in formalin-fixed, paraffin-embedded tissues: an enhancement method for immunohistochemical staining based on microwave oven heating of tissue sections. J Histochem Cytochem 39:741-8, 1991.

BioGenex

49026 Milmont Drive, Fremont, CA 94538, USA
Tel: +1 (800) 421-4149, Fax: +1 (510) 824-1490

e-mail: support@biogenex.com

780 Cailun Road, Suite 702, Zhangjiang Hi-Tech Park
Pudong New Area, Shanghai 201203, P.R. China
Tel: +86-21- 51371830, Fax: +86-21- 51371829

225 Cherlapally, Industrial Development Area, Phase II
Hyderabad, Andhra Pradesh, India 500051
Tel: +91-40-27185500, Fax: +91-40-27185511/21